

BIOLIS 24i Premium

Improved ease of operation and better test efficiency,
and provides the optimal solution for routine operations,
STAT operation and a broad range of clinical analysis requirements.



The newest Software makes it possible to take usability a step further

- The renewal of user friendly software**
New convenient functions are available, such as check for inappropriate calibration results by flag indication compared to the previous results, and capability of using multiple bottles of the same reagent item on the same tray, etc.
- Friendly user-interface**
Real-time display of analysis process, reagent residual volume and so on.
- Sharp reduction of reaction volume**
Improvements to the optical measurement system have reduced reaction volume by maximum 30% compared to the former model, BIOLIS 24i.
- Sample volume reduction**
Analyzer can test small volume samples (around 50 µl), effectively, such as those obtained from the elderly or infants.
- Night shift mode operation**
Night shift mode enables simpler test.
- Simple reagent management**
Automatic reagent residual volume calculation system allows display of real-time residual volume percent.
- LAN connection**
Can be connected to an electronic chart of hospital LAN. (Option)
- Low noise**
Operation noise is reduced by around 20% compared to BIOLIS 24i.

BIOLIS 24i Premium


Accuracy features

- Individual probes are used for 1st and 2nd reagents to avoid carry-over.
- Sample pre-dilution ratio and auto-rerun dilution ratio can be set in advance.
- Carry-over protection program for reaction cuvettes and probes.
- Detecting function for contaminated reaction cuvettes are available.
- Reagent cooling is available.
- Reaction cuvette washing with heated water and 2 kinds of washing solution (alkaline and acidic) and heated water washing for all probes.

Operation features

- Sample cup and primary tube (5, 7, 10 ml) can be used.
- Auto-gain function when lamp is replaced.
- Reaction waste is stored in a dedicated tank.
- Test result data are compatible with BIOLIS 24i.
- Sampling-stop function using END barcode is available.
- Data base control function is available.
- Improved error list is displayed (compared to BIOLIS 24i).
- Liquid level alarm for water reservoir and waste liquid reservoir is available.

Main unit arrangement




Front view




- Control PCB
- Reagent pump
- Reagent probe inside washing pump
- Sample pump
- Sample probe inside washing pump
- Probe outside washing pump
- Pump for ISE
- Reagent tray temperature controller

Air pressure mixing system


Our original system for mixing the sample and reagent using air pressure alone.



ISE module (OPTION)



[A] ISE is Direct Method
• Throughput is 400 tests/hour including ISE.
• ISE module is equipped with BIOLIS 24i Premium.
• It makes easy to replace the electrode.




[B] Consumables of the ISE module as follows
• Calibrator
• Cleaning solution
• Electrode (Na, K, Cl, Ref)

Friendly User-interface

Renewal of User-interface
Readily understandable screen and ease of operation


New functions
Automatic reagent residual volume calculation,
Current QC screen, Simple operation mode,
Soft sample tray mode, etc.




Run monitor screen




Order entry screen



Test results screen



QC screen



QC graph screen

Test items list

Clinical chemistry	LD(LDH)	AST(GOT)	ALT(GPT)	ALP	γ-GTP	CK(CPK)	LAP	ChE
	AMY	P-AMY	BUN	CRE	UA	T-CHO	TG	HDL-C
	LDL-C	TP	ALB	CRE	IP	Mg	Fe	UIBC
	D-BIL	T-BIL	GLU	NEFA	PL	ALD	SIA	TBA
Immuno-assay	GA	CK-MB	Apo A-I	Apo A-II	Apo B	Apo CII	Apo CIII	Apo E
	Lp(a)	RF	ASO	TPLA	RPR	β2-m	Tf	C3
TDM	CRP	IFG-G	lg-A	lg-M	lg-E	HAL	PB	
	C4	CBZ	CSA	DIG	DIGT			
Coagulation	PHT	THEO	VPA					
	ATIII	APL	PLG	PC	FDP	D dimer		
ISE	Na	Cl	K					

BIOLIS 24i Premium

Specification

Analysis	System type	Discrete single line random access multi-test analysis	
	Number of test items on board	36+3 (ISE) or 24+3 (ISE)	
	Throughput	240 tests/hour, 400 tests/hour including ISE	
	Analysis method	End point assay, rate assay, ISE (option)	
Calibration curve		9 kinds (linear, spline, etc)	
Sample	Sample kind	Serum, plasma, urine, CSF	
	Sample container	Sample cup, primary tube (5, 7, 10 ml)	
	Number of sample on board	Maximum 55/tray	
	Number of sample tray	Maximum 10	
	Sample dispense volume	2.0- 30.0 µl (0.1 µl step)	
	Dilution ratio	6, 10- 100	
	STAT sample	Available	
	Sample barcode	Barcode reader supplied as an option	
Reagent	Number of bottles on board	72 (36 items) or 48 (24 items)	
	Bottle volume	13, 25, 40 ml or 20, 40, 60 ml	
	Reagent dispense volume	20- 330 µl (1 µl step)	
	Reagent storage	Reagent cooling available	
	Residual volume	Level sensing or count down calculation	
	Reagent tray	36-sectors or 24-sectors (Removable)	
	Reagent barcode	Barcode reader supplied as an option	
	Reaction	Cuvette material	Plastics (Semi-disposable)
		Reaction volume	Minimum 140 µl, maximum 400 µl
		Reaction time	10 min. (1st reaction 5 min., 2nd reaction 5 min.)
Reaction temperature		37.0 ± 0.1°C	
Optical measurement		12 fixed wavelengths (340-800 nm)	
Optical source		Tungsten halogen lamp (Long-life type)	
Optical range		OD 0- 2.5	
Cuvette washing		Auto washing with heated water and 2 kinds of washing solution	
Reaction waste collection		Reaction waste to be stored in a dedicated tank	
Pure water consumption		Maximum 3.5 l/hour	
User interface	Run monitor	Analyzer operation status display	
	Reaction curve monitor	Optical absorbance graphic display	
	QC	Based on Westgard's algorithm, etc.	
	Voice message	Available	
	Test results storage	10,000 samples maximum	
	Printer	Internal, External (option)	
	Power supply	AC 100/115/230 volt (50/60 Hz)	
	Environmental conditions	Ambient temperature 15- 30°C, Humidity 30- 80% (No condensation)	
Dimension	800 mm (w) x 670 mm (d) x 520 mm (h)		
Weight	Approx. 95kg		
Option	ISE module		
	Sample barcode reader, Reagent barcode reader		
	Water purification system		
	External printer		