

# CLEAR, PURE AND FAST SAMPLES




Helping all people  
live healthy lives


**BD Vacutainer<sup>®</sup>**  
**SST<sup>™</sup> II Advance**

**BD Gel Technology**  
**Keeps on Moving**


**Safety**  
BD Hemogard™  
Safety Closure



## BD SST™ II Advance



**Laboratory Workflow & Process Efficiencies**  
Unique gel design



**Sample Integrity**  
New acrylic based gel

## Your Concern

## Our Solution

### SAFETY

#### Healthcare Worker

- Accidental exposure to blood can lead to blood-borne pathogen infection such as HIV, HBV, HCV

#### Patient

- Patient identification errors & incorrect test results can lead to inappropriate and hazardous treatment

#### Increased Healthcare Worker Safety

- Reduced exposure to blood through the virtually unbreakable medical grade PET tubes and the renowned BD Hemogard™ Safety Closure

#### Increased Patient Safety

- Proven tube performance through clinical trials ensuring optimal patient results
- No need to aliquot reducing patient identification errors

### SAMPLE INTEGRITY

#### Sample Purity & Analyte Stability

- Impure samples and analyte interference (direct contact with cells, ordinary gel, etc.) can lead to:
  - Analyzer interference
  - Re-testing
  - Extra handling
  - Delayed reporting of patient results & hence patient treatment

#### Enhanced Sample Purity & Analyte Stability

- No effect on serum AST, LD, K+ for 6 days at 4°C
- More than 90% recovery with the most commonly used therapeutic drugs up to 7 days at 4°C
- May be recommended for the most commonly tested hormones up to 7 days at 4°C
- May be frozen at -20°C under certain conditions

### LABORATORY WORKFLOW & PROCESS EFFICIENCY

#### Laboratory Workflow

- A slow turnaround time (TAT) can lead to delay in patient treatment

#### Laboratory Process

- Extra handling & labour can result in:
  - Slower TAT
  - Increased labour & material cost
  - Patient identification errors

#### Advanced Laboratory Workflow

- Reliable gel movement performance and barrier formation under stressed centrifugation conditions of 800xg at 10 minutes
- May be centrifuged along a g-force continuum between 1800g-3000g using a 4 or 5 minute spin<sup>1</sup>

#### Advanced Process Efficiency

- Primary tube sampling and storage without the need to aliquot

### MINIMISING PREANALYTICAL VARIABLES

#### Laboratory Errors in the Preanalytical Phase

- 50% of all laboratory errors occur in the preanalytical phase<sup>2</sup>
- 11% of patients are mismanaged<sup>2</sup>
- 15% of patients undergo further unnecessary treatment<sup>2</sup>

#### Reduction in Preanalytical Variability

- Providing a reliable barrier between the serum and the clot
- Offering a gel separator with minimum adsorption properties

1) 1300-2000xg for 10 minutes at room temperature remain the optimal centrifugation recommendations for general use.

2) M Plebani et al: Laboratory errors: their frequency, type and causes. The Impact of the Preanalytical Phase on the Quality of Laboratory Results. Oxford, July 1996.



# BD Gel Technology Keeps on Moving

To meet the needs of the changing clinical environment, BD offers the widest range of clinically documented sample separation solutions.

## The BD Vacutainer® SST™ II Advance through greater reliable barrier formation provides:

- Optimal Patient Results
- Enhanced Sample Integrity
- Advanced Laboratory Workflow & Process Efficiency
- Reduced Preanalytical Variability

	Catalogue Number	Draw (ml)	Tubing Size (mm)	Label
	366882	2.5	13 x 75	Block
	367957	3.5	13 x 75	Paper
	368498	3.5	13 x 75	See Through
	368879	4	13 x 100	See Through
	367955	5	13 x 100	Paper
	366566	5	13 x 100	See Through
	367953	8.5	16 x 100	Paper
	366644	8.5	16 x 100	See Through

## Reorder Information

Packaging: 100/box - 1000/case

All products available with **18-month shelf life.**

Different closure colours available upon request

BD Diagnostics - Preanalytical Systems  
 P.O. Box 52279, Dubai, United Arab Emirates  
 Tel : +971 4 3379525, Fax : +971 4 3379551  
 E-mail: EMA\_PAS@bd.com

www.bd.com

BD and BD Logo are trademarks of Becton, Dickinson and Company. ©2015





# BD Vacutainer® SST™ II Advance

## BD Gel Technology Keeps on Moving



### TDM Fact Sheet

BD Has selected these drugs as representatives from each class.

Drug Class	Time Temperature	4 hrs	8 hrs	24 hrs	48 hrs	7 days
		25°C	25°C	25°C	25°C	4°C
Cardioactive	Digoxin	•	•	•	•	•
	NAPA	•	-	•	•	•
Antibiotic	Procainamide	•	-	•	•	•
	Gentamicin	•	-	•	•	•
	Tobramycin	•	•	•	•	•
Antiepileptic	Vancomycin	•	•	•	•	•
	Carbamazepine	•	•	•	•	•
	Phenobarbital	•	•	•	•	•
Bronchodilator	Phenytoin	•	•	•	•	•
	Valproic Acid	•	•	•	•	•
	Theophylline	•	•	•	•	•
Analgesic	Acetaminophen	•	-	• <sup>1</sup>	•	•
	Salicylate	•	-	•	•	•

1. Acetaminophen was stable at 24 hrs, 25°C in one of two studies. Refer to BD White Papers for full details.

- Stable      - Not Tested

The table above indicates within-tube stability vs initial time ( $t_0$ ) in BD Vacutainer® SST™ II Plus Tubes. Refer to BD White Papers for full details.

Ref: BD White Paper VS5776, VS7050.

### The absorption of drugs by gel is dependent upon several factors<sup>2</sup> including:

- The chemical nature of the gel
  - The chemical nature of the drug
  - Specimen volume on gel
  - Time specimen sits on the gel
  - Storage temperature
2. Dasgupta, et al. "Absorption of Therapeutic Drugs by Barrier Gels in Serum Separator Blood Collection Tubes". Am J Clin Pathol. 1994; 101: 456-461.

### Special Chemistry Fact Sheet

BD has selected these analytes as representatives from each type of compound.

	Time Temperature	8 hrs	24 hrs	48 hrs	7 days
		25°C	25°C	4°C	4°C
Proteins/Peptides	CEA	•	-	•	•
	CRP	•	-	•	•
	Ferritin	•	-	•	•
	FSH	• <sup>1</sup>	-	•	• <sup>1</sup>
	LH	•	-	•	•
	Maternal AFP	• <sup>1</sup>	-	•	•
	PSA	-	• <sup>2</sup>	-	-
	Free T4	• <sup>1</sup>	• <sup>1</sup>	-	• <sup>1</sup>
	Total T3	•	-	•	•
	Total T4	•	-	•	• <sup>1</sup>
Vitamins	TSH	•	•	•	•
	Folate	• <sup>1,3</sup>	• <sup>1,3</sup>	• <sup>1,3</sup>	• <sup>1,3</sup>
Steroids	Vitamin B12	•	•	• <sup>1</sup>	• <sup>1</sup>
	Cortisol	•	•	-	•
	Estradiol	-	•	•	•
	Progesterone	-	•	•	•
	Testosterone	•	• <sup>1</sup>	-	• <sup>1</sup>

1. The mean bias vs initial time ( $t_0$ ) was within the clinical acceptance limit but the 95% limit for the mean bias was not. Refer to BD White Papers for full details.

2. PSA was stable at 24 hrs, 4°C

3. Stability of folate may vary depending on storage conditions.

- Stable      - Not Tested

The table above indicates within-tube stability vs initial time ( $t_0$ ) in BD Vacutainer® SST™ II Plus Tubes, except for progesterone and estradiol. Data for progesterone and estradiol indicates within-tube stability vs  $t_0$  in BD Vacutainer® SST™ II Advance Tubes. Refer to BD White Papers for full details.

Ref: BD White Paper VS5778, VS7051, VS7192

Whenever changing any manufacturer's blood collection tube type, size or storage condition for a particular laboratory assay, the laboratory personnel should review the tube manufacturer's data and their own data to establish/verify the reference range for a specific instrument/reagent system. Based on such information, the laboratory can then decide if a change is appropriate.

