

Mission® Urinalysis Reagent Strips and Urine Analyzers



Obtain reliable and cost-effective results with *Mission*® Urinalysis Reagent Strips and Urine Analyzers!

- Accurate
- Reliable
- Convenient

Global Diagnostics for Local Markets™

Urinalysis Reagent Strips



Simple and Accurate


- Analytical sensitivity better than or comparable to market leaders
- High quality color chart ensures accurate visual reading

Flexible


- Compatible for visual and analyzer reading
- More than 30 different combinations available

Multiple Packaging Options and Long Shelf Life


- Available in 25, 50, 100 and 150 strips per canister
- 2 year shelf life for unopened canisters which offers cost savings and convenience for high volume testing
- 3 month shelf life for strips in opened canisters



Step 1: Immerse strip into urine



Step 2: Remove excess urine



Step 3: Obtain results by analyzer or visual reading

Catalog No.	No. of Parameters	Strips per Canister [◇]	Pouch Packaging [▲]	Visual-Read and Analyzer-Read Strips	Analyzer-Read Strips: Standard (S) or Additional (A)	Parameters																		
						ASC	GLU	BIL	KET	SG	BLO	pH	PRO	URO	NIT	LEU								
U031-111	11	100	✓	✓	S	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*			
U031-101	10	100	✓	✓	S	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*			
				✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
U031-091	9	100	✓	✓	S	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*			
U031-081	8	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*			
				✓	S	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
				✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
U031-071	7	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*				
U031-061	6	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*			
				✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
U031-051	5	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
U031-041	4	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
U031-031	3	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
U031-021	2	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
U031-011	1	100	✓	✓	A	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*		
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
				✓		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*

Visual/U120 Strip Size
1-6 Parameters: 5 mm x 80 mm
7-11 Parameters: 5 mm x 110 mm

U500 Strip Size
1-11 Parameters: 5 mm x 110 mm

◇Also available in canisters of 25, 50 and 150 strips

Cleared for US 510(k); CLIA waived
CE Marked for sale in the European Community

