



CE-Marked Blood Grouping Reagents

Our range of CE-Marked Blood-Grouping reagents are available from 2 clones , ideal for markets that require testing from 2 different formulations. Choose which clone is best for you using our handy table below.



- ✓ CE-Marked
- ✓ Manufactured in the UK
- ✓ Rapid Delivery at Low Prices



CE-Marked Blood Grouping Range

Blood Type / Product	BG-Reagent	BC-Reagent
Anti-A	Clone BIRMA-1	Clone 9113D10
Anti-B	Clone LB-2	Clone 9621A8
Anti-A,B	Clones ES-4 & ES-15	Clones 152D12, 9113D10 & ES15
Anti-D	Clones RUM-1 & MS-26	Clones RUM-1 & MS-26
Anti-A,B,D Kit	Clones listed above	Clones listed above
Anti-A,B,AB, D Kit	Clones listed above	Clones listed above
Anti-AHG	-	Clones BRIC-8 & BIRC

Which CE-Marked Blood Grouping Reagents are best for me?

	BG- Reagent	BC- Reagent
CE-Marked	Yes	Yes
Lead Time	Immediate Despatch	6-8 Weeks
Lower Price	Yes	No
Suitable Methods	Tube, Diamed-ID & Slide Techniques	Tube, Diamed-ID, Ortho BioVue, Microplate & Slide Techniques

Rapid Labs Limited



Address: Unit 2 & 2A, Hall Farm Business Centre ,
Church Road, Little Bentley, Colchester,
Essex CO7 8SD, United Kingdom



+ 44 01206 250484



medical@rapidlabs.co.uk



www.rapidlabs.co.uk